

Recommended Supplies for Lammina Veneer Treatment



The following are supplies needed for bonding veneers using Lammina trays.

Material	Amount/Size	Brand	Source
Isolation	small or regular	OptraGate	https://www.ivoclar.com/en_us/shop/p/clinicalaccessories/optragateregularrefill80/p/b601513
Prophy Paste	medium, fluoride-free	Ivoclar Proxyt®	https://www.ivoclar.com/en_us/shop/p/prevention/proxytmediumfluoridefree1x80g/p/701472
Etching Gel	37% phosphoric acid	Total Etch	https://www.ivoclar.com/en_us/shop/p/adhesives/totaletchrefill2x2g/p/550588an
Isolation Tape	0.10 or 0.07 thickness 18mm width	TDV IsoTape	https://www.net32.com/ec/tdv-isotape-5-m-tape-isolation-indications-d-144539
Adhesive		3M Scotchbond Universal Plus	https://www.3m.com/3M/en_US/dental-us/featured-products/scotchbond-universal-plus/
Flowable Composite		G-aenial Universal Injectable CLEARFIL Majesty ES Flow	https://www.gc.dental/america/products/operator/composite-restoratives/g-aenial-universal-injectable https://kuraraydental.com/product/clearfil-majesty-es-flow/
Glycerin Gel		Liquid Strip SDS	https://www.ivoclar.com/en_us/shop/p/cementation/liquidstriprefill1x25g/p/532505an
Interproximal Strips	extra-fine (yellow)	Komet USA	https://www.kometusa.com/Products/General-Dentistry/Category/Diamonds/Diamond-Strips/DS25EF
Resin Color Modifier		Kolor + Plus™	https://www.kerrdental.com/kerr-restoratives/kolor-plus-resin-color-modifier-and-opaquing-kit
Scalpel	#12		

Additional Supplies: The following are some suggested supplies you may need for finishing & polishing.

- Diamond Bur (Pointed) - fine grit (red) & extra fine grit (yellow)
- Finishing Carbide Bur - flame-shaped
- Diamond Bur (Oval-shaped) - fine grit (red)
- Polishing Discs - extra-fine grit
- Spiral Polisher
- Aluminum Oxide Paste - 3 micron & 1 micron
- Polishing Brush - cup-shaped, goat-hair
- Polishing Brush (Cotton)
- High-Gloss Polishing Paste
- Polishing Bur